



WOMEN'S HEALTH

INTRODUCTION	2
OVUM PICKUP NEEDLES	
ACE Single Lumen Ovum Pickup Needle	3
ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	4
ACE-M Manual Single Lumen Ovum Pickup Needle	5
OVUM PICK Double Lumen Ovum Pickup Needle	6
EMBRYO TRANSFER CATHETERS	
EMBRYO TRANS Embryo Transfer Catheter	7
ECHO TRANS Embryo Transfer Catheter	8
TWINKLE EMBRYO TRANS Embryo Transfer Catheter	9
BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip	10
BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip	11
PRO EMBRYO TRANS Embryo Transfer Catheter	12
PRO ECHO TRANS Embryo Transfer Catheter	13
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter	14
INTRA UTERINE INSEMINATION CATHETERS	
INSEMINA Intra Uterine Insemination Catheter	15
IUI CURVED Intra Uterine Insemination Catheter Curved	16
TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	17
FLEXI Intra Uterine Insemination Catheter	18
STUMPY Intra Uterine Insemination Cannula	19
FALLOPIAN TUBE RECANALISATION CATHETERS	
SALPINX CATH Fallopian Tube Recanalisation Catheter	20
HYSTEOSALPINGOGRAPHY CATHETER	
HYSTERO CATH Hystero Salpingography Catheter	21
ENDOMETRIAL BIOPSY SAMPLING SYSTEM	
THE SAMPLER Endometrial Biopsy Catheter	22
PRO-VAC Endometrial Biopsy Sampler	23
VITRIFICATION STRAWS	
SURESTRIP Vitrication Straws (Open System Vitrication)	24
SURELOCK Vitrication Straws	25
SUREHOOK Vitrication Straws	26
OTHER PRODUCTS	
ASPI Aspiration Needle	27

Product Warranty

Allwin warrants that at the time of manufacture, the products are prepared and tested in accordance with highest standards. Due to biological differences in individuals, no product is 100% effective under all circumstances for all patients. Since we have no control of the condition under which the product is being used, the diagnosis, the method, the administration, the handling, and the storage, we do not warrant any consequences which may arise pertaining to the final use of the device. The foregoing warranty is exclusive and in lieu of all other warranties written, oral or implied. No representative of **Allwin** may change any of the foregoing and the buyer accepts the product subject to all terms hereof.

Design Changes

Continuous improvement is the chosen path at **Allwin** and we reserve the right to change or discontinue any product without notice.

Return Policy

Allwin will accept returned goods only after issuing a Return Goods Authorisation (RGA). The return request must provide a reason for return, invoice number and order number. Credit will be issued once the request is evaluated as per company rules. We will accept a return of the stock item within 30 days of the original invoice date. We will not accept the return of a nonstock, opened, and damaged products. Returned stock items may be subject to a restocking fee.


Product testing

Allwin products are endotoxin tested with a limulus amebocyte lysate (LAL) assay unless otherwise noted. Our Ovum Pickup needles and Embryo Transfer catheters are batch tested for toxicity with a one-cell mouse embryo assay (MEA). And our Intra Uterine Insemination catheters are Human Sperm Survival Assay (HSSA) tested.

Storage of Products


Allwin recommends that products be stored in a cool and dry place. Plastic material used may degrade following extended periods of exposure to light and temperature.

ACE Single Lumen Ovum Pickup Needle




Allwin ovum pickup needles are used for laparoscopic or ultrasound guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line




Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35	ACE Single Lumen Ovum Pickup Needle	35	16	10
A-OVM-SL17/35	ACE Single Lumen Ovum Pickup Needle	35	17	10
A-OVM-SL18/35	ACE Single Lumen Ovum Pickup Needle	35	18	10

*All products may not be available in all countries. Please contact our sales department for details.




Highly echogenic tip

ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line




Allwin ovum pickup needles are used for laparoscopic or ultrasound guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line
- Alternate aspiration & flushing made possible through single channel
- Multiple operation functionality by using stopcock valve




Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	16	10
A-OVM-SL17/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	17	10
A-OVM-SL18/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	18	10

*All products may not be available in all countries. Please contact our sales department for details.



Highly echogenic tip

ACE-M Manual Single Lumen Ovum Pickup Needle



Allwin ovum pickup needles are used for laparoscopic or ultrasound guided transvaginal retrieval of oocytes from ovarian follicles.

- When used with ultrasonic equipment, echo tip enhances visualization.
- Designed handle for greater control during ovum retrieval.
- Luer Lock connector for connecting a syringe
- Amount of suction applied is manually controlled.


Tri-facet geometry
with highly echogenic tip



Ergonomic
handle


Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	35	16	10
A-OVM-SL17/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	35	17	10
A-OVM-SL18/35-M	ACE-M Manual Single Lumen Ovum Pickup Needle	35	18	10

*All products may not be available in all countries. Please contact our sales department for details.




Highly echogenic tip

OVUMPICK Double Lumen Ovum Pickup Needle



Allwin ovum pickup needles are used for laparoscopic or ultrasound guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for tight fit
- Non-coring needle
- With luer lock connector
- Outer lumen allows to flush follicles



Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
OVM-DL16/35	OVUMPICK Double Lumen Ovum Pickup Needle	35	16	10
OVM-DL17/35	OVUMPICK Double Lumen Ovum Pickup Needle	35	17	10

*All products may not be available in all countries. Please contact our sales department for details.




Highly echogenic tip

EMBRYO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Soft inner catheter with smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately




Order No.	Description	Catheter Length	Qty. / Box
EB-ETC23	EMBRYO TRANS Embryo Transfer Catheter	23	10
EB-ETC18	EMBRYO TRANS Embryo Transfer Catheter	18	10
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10

*All products may not be available in all countries. Please contact our sales department for details.




ECHO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.


- Soft inner catheter with smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Visible under ultrasound
- Stylet for outer sheath available separately



Stylet for Embryo Trans / Echo Trans Catheter / Twinkle Embryo Trans

Order No.	Description	Catheter Length	Qty. / Box
ET-ETC23	ECHO TRANS Embryo Transfer Catheter	23	10
ET-ETC18	ECHO TRANS Embryo Transfer Catheter	18	10
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10

*All products may not be available in all countries. Please contact our sales department for details.




TWINKLE EMBRYO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Soft inner catheter
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately



Stylet for Embryo Trans / Echo Trans Catheter / Twinkle Embryo Trans

Order No.	Description	Catheter Length	Qty. / Box
TEB-ETC23	TWINKLE EMBRYO TRANS Embryo Transfer Catheter	23	10
TEB-ETC18	TWINKLE EMBRYO TRANS Embryo Transfer Catheter	18	10
STY23	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18	Stylet for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10
STY23-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	23	10
STY18-OS	Stylet with outer sheath for Embryo Trans / Echo Trans / Twinkle Embryo Trans Catheter	18	10

*All products may not be available in all countries. Please contact our sales department for details.




BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- Soft inner catheter with smooth tip
- Marking on inner steel pipe and outer sheath
- Bulb shape tip of outer sheath
- Silicone depth positioner for precise placement




Styler for BULB TRANS Embryo Transfer Catheter-Bulb Tip

Order No.	Description	Qty. / Box
BT-ETC-U	BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip	10
STY-BU	Styler for BULB TRANS Embryo Transfer Catheter - Bulb Tip	10

*All products may not be available in all countries. Please contact our sales department for details.




BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- Soft inner catheter
- Marking on inner steel pipe and outer sheath
- Bulb shape tip of outer sheath
- Silicone depth positioner for precise placement
- The catheter is curved to facilitate insertion
- Echogenic tip of inner catheter to improve the visibility under ultrasound



Order No.	Description	Qty. / Box
BT-ETC-U-S	BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip	10
STY-BU	Stylet for BULB TRANS Embryo Transfer Catheter - Bulb Tip	10

*All products may not be available in all countries. Please contact our sales department for details.




PRO EMBRYO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible transfer catheter
- Silicone depth positioner for precise placement
- Orientation arrow on hub



Soft obturator for Pro Embryo Trans / Pro Echo Trans Catheter /Twinkle Pro Embryo Trans Catheter



Stylet for Pro Embryo Transfer / Pro Echo Trans Catheter /Twinkle Pro Embryo Trans Catheter



Order No.	Description	Qty. / Box
PEB-ETC	PRO EMBRYO TRANS Embryo Transfer Catheter	10
PEB-ETC-SB	PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
PEB-ETC-STY	PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Transfer / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.




PRO ECHO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

Allwin embryo transfer catheters (pro embryo trans, pro echo trans, and twinkle pro echo trans) also include stylets or soft obturators intended to assist uterine access of the guide catheter during an embryo transfer procedure.

- Soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible transfer catheter
- Visible under ultrasound
- Orientation arrow on hub



Soft obturator for Pro Embryo Trans / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans

Stylet for Pro Embryo Transfer / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans

Order No.	Description	Qty. / Box
PET-ETC	PRO ECHO TRANS Embryo Transfer Catheter	10
PET-ETC-SB	PRO ECHO TRANS Embryo Transfer Catheter with Soft Obturator	10
PET-ETC-STY	PRO ECHO TRANS Embryo Transfer Catheter with Stylet	10
SB-PET-ETC	Soft obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Transfer / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.



TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter



Allwin embryo transfer catheters are used to introduce in vitro fertilized (IVF) embryos into the uterine cavity.

- Soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible transfer catheter
- Silicone depth positioner for precise placement
- Orientation arrow on hub



Soft obturator for Pro Embryo Trans / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans Catheter



Stylet for Pro Embryo Trans / Pro Echo Trans Catheter / Twinkle Pro Embryo Trans Catheter




Order No.	Description	Qty. / Box
TPEB-ETC	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter	10
TPEB-ETC-SB	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
TPEB-ETC-STY	TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Transfer / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.




INSEMINA Intra Uterine Insemination Catheter



Allwin's intra uterine insemination catheters are intended for the introduction of washed spermatozoa into the uterine cavity.

- Soft uterine catheter with firm support
- Integrally moulded hub
- Silicone depth positioner for precise placement
- Open tip and closed tip with two lateral holes options available



Order No.	Description	Qty. / Box
I-IUI	INSEMINA Intra Uterine Insemination Catheter	25
I-IUI-O	INSEMINA Intra Uterine Insemination Catheter (Open Tip)	25

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip




Open Tip

IUI CURVED Intra Uterine Insemination Catheter Curved



Allwin's intra uterine insemination catheters are intended for the introduction of washed spermatozoa into the uterine cavity.

- Soft uterine catheter with firm support
- Integrally moulded hub
- Silicone depth positioner for precise placement
- Curve shaped catheter
- Open tip and closed tip with two lateral holes options available



Order No.	Description	Qty. / Box
I-IUI-C	CURVED Intra Uterine Insemination Catheter	25
I-IUI-C-O	CURVED Intra Uterine Insemination Catheter (Open Tip)	25

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip




Open Tip

TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet



Allwin's intra uterine insemination catheters are intended for the introduction of washed spermatozoa into the uterine cavity.

- Soft uterine catheter with firm support
- Integrally moulded hub
- Silicone depth positioner for precise placement
- Malleable stylet to aid insertion of the catheter
- Open tip and closed tip with two lateral holes options available



Order No.	Description	Qty. / Box
T-IUI	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	25
T-IUI-O	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip)	25

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip




Open Tip

FLEXI Intra Uterine Insemination Catheter



Allwin's intra uterine insemination catheters are intended for the introduction of washed spermatozoa into the uterine cavity.

- Shapeable outer sheath for easy insertion
- Centimeter marking on catheter for precise placement
- Stylet available separately
- Open tip and closed tip with two lateral holes options available



Order No.	Description	Qty. / Box
F-IUI	FLEXI Intra Uterine Insemination Catheter	25
F-IUI-O	FLEXI Intra Uterine Insemination Catheter (Open Tip)	25
STY18	Stylet for FLEXI IUI Intra Uterine Insemination Catheter	25

*All products may not be available in all countries. Please contact our sales department for details.



Close Tip



Open Tip

STUMPY Intra Uterine Insemination Cannula



Allwin's intra uterine insemination catheters are intended for the introduction of washed spermatozoa into the uterine cavity.

- Conical proximal end
- Luer lock connector
- Highly flexible cannula that
- Small shaped catheter




Order No.	Description	Qty. / Box
ST-IUI	STUMPY IUI Intra Uterine Insemination Cannula	25

*All products may not be available in all countries. Please contact our sales department for details.




SALPINX CATH Fallopian Tube Recanalisation Catheter



Used for fallopian tube recanalization under radiological fluoroscopy control.


- Preshaped 5 fr catheter
- Centimeter marking on catheter for precise placement
- Malleable stylet
- Hydrophilic guidewire available separately



Order No.	Description	Qty. / Box
SC-FTR-M	SALPINX CATH Fallopian Tube Recanalisation Catheter with Malleable Stylet	5
SC-FTR-MG	SALPINX CATH Fallopian Tube Recanalisation Catheter with Malleable Stylet and Hydrophilic Guidewire	1


*All products may not be available in all countries. Please contact our sales department for details.

HYSTERO CATH Hysterosalpingography Catheter



Used for injection of contrast media into the uterine cavity and fallopian tube(s) to evaluate patency of the fallopian tubes.

- Strong low profile 7 fr balloon
- Closed end with side hole and open end HSG catheter available
- Stopcock design for locking of balloon channel and regulate flow of contrast media.
- Stylet provided for added rigidity during innodictopm
- 3 ml syringe included



Order No.	Description	Qty. / Box
HC-HSG7	HYSTERO CATH Hysterosalpingography Catheter	5

*All products may not be available in all countries. Please contact our sales department for details.

Caution:

This product contains natural rubber latex which may cause allergic reaction in some individuals.

THE SAMPLER Endometrial Biopsy Catheter



Used to obtain histological tissue biopsy and extract sample of uterine menstrual content. Used for obtaining endometrial tissue.

- Centimeter marking provided on the catheter for precise placement
- Smooth rounded tip

Smooth
rounded
tip


Centimeter
marking



Order No.	Description	Qty. / Box
SEBC	THE SAMPLER Endometrial Biopsy Catheter	25


*All products may not be available in all countries. Please contact our sales department for details.

PRO-VAC Endometrial Biopsy Sampler



Used for obtain histological tissue biopsy and extract sample of uterine menstrual content. Used for obtaining endometrial tissue.

- Centimeter marking provided on the catheter for precise placement
- Smooth rounded curved tip
- Locking syringe
- Angled tip



Order No.	Description	Qty. / Box
PVES	PRO-VAC Endometrial Biopsy Sampler	25


*All products may not be available in all countries. Please contact our sales department for details.

SURESTRIP Vitrification Straws (Open System Vitrification)



Used for obtain histological tissue biopsy and extract sample of uterine menstrual content.

- High survival and development rates.
- Secure coupling.
- Marking on tip extreme for easy visibility under microscope.
- Adequate flat surface for writing.
- Tip width ensuring placement of minimum solution
- 1-cell MEA $\geq 80\%$ expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours



Colour	White	Red	Ivory	Blue	Yellow	Qty. / Box
Order No.	SUS-W	SUS-R	SUS-I	SUS-B	SUS-Y	10


*All products may not be available in all countries. Please contact our sales department for details.

SURELOCK Vitrification Straws



Used for obtain histological tissue biopsy and extract sample of uterine menstrual content.

- Secure air tight coupling.
- Semi transparent body and cap.
- Arrow on top surface ensuring orientation during use.
- Shallow well on tip distal end for exact placement of embryo.
- Ease in visibility under microscope
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of in 1-cell stage of embryos
- Available in 5 colours



Colour	Pink	Yellow	Green	Blue	Violet	Qty. / Box
Order No.	SUL-P	SUL-Y	SUL-G	SUL-B	SUL-V	10


*All products may not be available in all countries. Please contact our sales department for details.

SUREHOOK Vitrification Straws



Used for obtain histological tissue biopsy and extract sample of uterine menstrual content.

- Hook shaped tip design
- Secure coupling
- Adequate flat surface for writing
- Hook width ensuring minimum placement of solution
- Flexible tip
- Sterility: SAL 10⁻⁶ - Gamma sterilized
- Vitrification of oocytes & embryos
- Available in 5 colours



Colour	White	Red	Green	Blue	Yellow	Qty. / Box
Order No.	SUH-W	SUH-R	SUH-G	SUH-B	SUH-Y	10


*All products may not be available in all countries. Please contact our sales department for details.

ASPI Aspiration Needle




Used for laparoscopic or ultrasound guided transvaginal aspiration.

- Bevel tip
- Centimeter marking provided for precise placement
- Transparent luer lock connector



Order No.	Description	Size cm	Gauge	Qty. / Box
ASPI 18/30	ASPI Aspiration Needle (30 cm) (18 Gauge)	30	18	10
ASPI 22/30	ASPI Aspiration Needle (30 cm) (22 Gauge)	30	22	10



Female Reproductive System

[illegible]

allwin[®]

Medical Devices



US Sales

✉ sales.usa@allwinmedical.com

☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com

☎ +91 9987084563

Worldwide sales (Excluding India/US)

✉ sales@allwinmedical.com

☎ +91 95949-81581